



# Progressive Compound Management Strategies and Pitfalls

Rhett Affleck, PhD  
LRIG  
June, 2009

# NEXUS Technology Platforms

## Combinatorial Chemistry Automation

**IRORI™**


RF MicroKan Synthesis Products  
2D NanoKan & X-Kan Synthesis Products →



## Protein Crystallography

**Crystal Farm™**


← Automated Crystal Incubation & Imaging



## Sample Management Solutions


**Universal Store™**

← -20C, Nitrogen  
Automated High-Speed Cherry Picking & Defrosting  
Decapping and Desealing solutions



**BioStore**

-80°C  
→ Automated picking and retrieval



# Best Storage Conditions

- **Ideal Storage for Stability**
  - Cold, inert, dry
  - Solid phase
    - Dry film, powder
- **Compound transfer is difficult**
  - Weighing
    - Slow
    - Doesn't work for dry films
    - Impractical for screening amounts
      - (1 – 1000 ng)
  - Volatile solvent transfer
    - Dissolve, transfer, evaporate
    - Slow
  - Apply to long term storage of compound reserve
  - Create a short term **Working Store**
    - (compounds in solution)


# Working Store

## Compound Storage and Transfer

- **Transfer solution-phase compound**
  - **Store dry:** Volatile solvent transfer
    - Store dry films
    - Dissolve – transfer - dry
  - **Store solutions:** Nonvolatile solvent transfer and storage
    - Flexible, fast, easy
- **Practical to store and transfer in DMSO**
  - DMSO solutions are a compromise
    - Some added reactivity
    - Complicated by its hygroscopic nature
      - Adds reactivity
      - Decreases compound solubility
  - For solution-phase storage is anything better than DMSO?
    - Universal solvent, plastic friendly, nonvolatile, decent biological compatibility
  - *Keep the water out!*

# Compound Concentrations

## Compound Storage and Transfer



# DMSO/Water

- **Freezing point**
  - 0 °C at 100% water
  - 18 °C at 100% DMSO
  - -73 °C at 33 mol % DMSO
    - 33% v/v water
  - -20 °C at 55 mol % DMSO
    - 17% v/v water




Fig. 2. Phase diagram for the water-DMSO system (plotted on a mole per cent basis). Symbols: L, liquid phase; A,  $\text{H}_2\text{O}(\text{s})$ ; B,  $\text{DMSO}\cdot 3\text{H}_2\text{O}(\text{s})$ ; C,  $\text{DMSO}\cdot 2\text{-}1/2\text{H}_2\text{O}(\text{s})^*$ ; D,  $\text{DMSO}\cdot 1/2\text{H}_2\text{O}(\text{s})^*$ ; E,  $\text{DMSO}(\text{s})$ .  
\*Tentative assignment.

\*Rasmussen et al, Nature vol.220.1968, 1315-1317

# DMSO/Water

- Non-ideal behavior
- 1DMSO:2water complexes
  - 33% v/v water

Characterization of Binary Solvent Mixtures

*J. Org. Chem.*, Vol. 66, No. 17, 2001 5847





Figure 1. Plots of (a) surface tension,  $\gamma$ , (b) viscosity,  $\eta$ , (c) excess enthalpy of mixing,  $H_m^E$ , and (d) deviations from ideality for the chemical shifts of protons with respect to DMSO,  $\Delta\delta(\text{DMSO} - \text{H}_2\text{O})$ , vs the mole fraction,  $X$ , for the DMSO/water mixture.


\*Catalan et al, *J. Org. Chem.* 2001, 66, 5846-5852

# Compound Solubility in DMSO/Water






# Compound Solubility in DMSO/Water




# Compound Solubility in DMSO/Water

## Saturated Azodipyridine




# Compound Solubility in DMSO/Water


## Saturated Rhodanine




# Compound Solubility in DMSO/Water



# Compound Solubility in DMSO/Water



# Compound Solubility in DMSO/Water



# Working Store

## Compound Concentration (in DMSO)

- **In practice: 1 – 30 mM**
  - **Lower is better (for compound management)**
    - Solubility
      - 16% not soluble at 10 mM
      - 65,500 member drug-like compound set\*
    - Liquid handling easier ( $\mu\text{L}$  instead of nL)
  - **Higher is better (for biology)**
    - Minimize DMSO at screening concentrations (1  $\mu\text{M}$  – 10  $\mu\text{M}$ )
  - **Compromise at ~ 1 mM**
    - 1 mM  $\rightarrow$  10  $\mu\text{M}$ , 1% DMSO, or
    - $\rightarrow$  1  $\mu\text{M}$ , 0.1% DMSO, or
    - $\rightarrow$  50  $\mu\text{M}$ , 5% DMSO
- $\rightarrow$  If needed:  
evaporate DMSO**

\*Balakin et al, J. Biol. Screen. 9(1); 2004

# Working Store

## Freeze/Thaw Cycles


- **Freezing and thawing does not degrade compounds**
- **Freezing DMSO solutions can cause some compounds to come out of solution**
  - Equilibrium solubility at 25 °C does not change
  - Dissolution kinetics can change
    - More stable crystals can form
- **Compounds will re-dissolve!**
  - Make sure that thawed solutions are mixed well




# Intermediate Transfer Plates

- **Storage at remote sites**
  - Are remote sites equipped with storage, liquid handling, and personnel to take multiple aliquots?
    - Example: seal and unseal samples at point of liquid handling
- **Intermediate dilution required**
  - 1  $\mu\text{L}$  @ 20 mM in DMSO delivered
  - 10  $\mu\text{M}$  in buffer desired
  - Need to dilute 1:2000
  - For example:
    - Dilute 1:10 with 10  $\mu\text{L}$  buffer,
    - then transfer 1  $\mu\text{L}$  to assay plate and dilute to 200  $\mu\text{L}$


# Compound Solubility in DMSO/Water




# Compound Solubility in DMSO/Water




# Compound Solubility in DMSO/Water



# Compound Solubility in DMSO/Water



# Compound Solubility in DMSO/Water



# Intermediate Transfer Plates

- ~~Storage at remote sites~~
- → Compound handling only at compound management sites
- ~~Intermediate dilution required~~
- → Make Assay-Ready Plates

# Screening Concentration

- **In practice: 1 – 30  $\mu\text{M}$**
- **Solubility**
  - Amphora Discovery\*
    - ~85% of compounds below expected 50  $\mu\text{M}$
    - ~55% of compounds below 10  $\mu\text{M}$
    - Buffer with 5% DMSO
    - 2,000 member diversity library
  - Pfizer\*\*
    - ~95% of compounds below expected 200  $\mu\text{M}$
    - ~45% of compounds 25  $\mu\text{M}$  – 100  $\mu\text{M}$
    - ~50% of compounds below 25  $\mu\text{M}$
    - Water with 5% DMSO
    - 1,000 member diverse set

\*Popa-Burke et al, Anal. Chem. 76(24); 2004

\*\*Nakayama, San Diego HTS Discussion Group, Dec. 7, 2005



# Screening Concentration Too High?

- **Many compounds not soluble at 10  $\mu$ M!**
  - Best would be to measure concentration
    - Not practical for HTS
  - Anticipate insoluble compound
- **Insoluble and soluble, high-concentration dependent aggregates act as promiscuous inhibitors\***
  - Study used 30  $\mu$ M and 5  $\mu$ M solutions diluted from 10mM DMSO stock


\*B.Y. Feng et al, Nature ChemBio. 1(3), 2005

# Screening Concentration Too High?

**Why does this matter to compound management?**

- **Screening “demands” push DMSO stock concentration**
  - Lower screening concentrations could lower DMSO stock concentrations

# Compound Solubility in DMSO/Water



# Remove DMSO before Screen?

- **Deliver compounds in DMSO or dry?**
  - In DMSO
    - Advantage:
      - Rapid mixing
    - Disadvantage
      - Crashes out of solution unpredictably
      - Forms promiscuous aggregate activity?
  - Dry
    - Disadvantage
      - Slow dissolution, so good mixing needed
    - Advantages:
      - Assay can dictate % DMSO (if any)
      - Insoluble compound stays on well surface
      - Easy to transport (dry)
      - Working Store concentration can be low!

← Better for biology

← Better for compound management

# My Experience with Dry Films

- **Pharmacopeia**
  - Single bead libraries
    - Extracting compounds from beads
    - Assays run from dry films
  - 1536-well reformatting
    - 1  $\mu$ L assays
- **Discovery Partners**
  - $\mu$ ARCS
    - Reformat compounds onto 10,000 spot cards
    - Assays run as compounds dissolve from dry film into agarose gels


# Microplate Mixing

- **Sonication**
  - Variable effectiveness
- **SonicMan, Matrical**
  - Pros: high throughput
  - Cons: exposed samples, wash probes, heat
- **Advalytix**
  - 100MHz, 20  $\mu$ m,
  - Pros: Whole plate, no water bath
  - Cons: low energy surface acoustic waves, mixing near surface, needs hard plate
- **AFA, Covaris**
  - Pros: sealed samples, isothermal, non-contact, high energy
  - Cons: water bath




# Orbital Table Mixing

## 150 cycles/min



# Covaris AFA Treatment





# Summary

- **Keep water out of your DMSO samples**
  - Working store
  - Intermediate plates
    - Intentional or not
- **Store and manipulate compounds only at expert facilities**
  - Centralize
  - Properly equip and train remote sites
- **Consider:**
  - JIT, individual-use, assay-ready plates
  - Lowering concentrations of DMSO stock
  - Delivering compounds dry